LASAK

PORESORB-TCP

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area, radiotherapy in the operated area, pathological changes in the jaws or buccal mucosa in dental indications.

b) **Temporary:** Acute febrile disease, pregnancy, treatment with corti- costeroids or other drugs affecting calcium metabolism, cytostatic or anticoagulant treatment, extreme physical or mental stress.

c) **General medicinal:** Endocrinopathy and other diseases affecting calcium metabolism, uncontrolled diabetes mellitus, hematological diseases and bleeding disorders, systemic bone diseases, systemic rheumatic diseases, congenital and acquired immune disorders, cy- tostatic, immunosuppressive or anticoagulant treatment.

Undesirable effects: Not observed to date.

Interactions: Not observed to date.

Means of application: Prior to the application of PORESORB-TCP, the granules may be mixed with the patient’s own blood from the operation wound or withdrawn venous blood, with the patient’s own platelet rich plasma (PRP). Application and final shaping of the mixture at the site of the defect is possible using a suitable metallic instrument. The adapted and condensed ceramic is enclosed by the soft tissues using a hermetic suture. It is not advisable to saturate the granules with aqueous solutions (such as normal saline).

Operation technique: PORESORB-TCP should be applied by or under the supervision of a physician with surgical and/or periodon- tal training and experience in the use of biomaterials. The use of the PORESORB-TCP material requires gentle operation techniques

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and prevention of contamination of the material and the operation wound (for example by saliva in dentoalveolar surgery). Undisturbed healing of the PORESORB-TCP material requires close contact with the vital bone tissue and assurance of stability in the bone socket. Micromovements of the solid particles cause proliferation of con- nective tissue between the solid particles and their fibrous encapsu- lation. The opening of the soft tissues should be extensive, enabling good access to the operation field. It is necessary to uncover the vital bleeding bone tissue for subsequent healing with the ceramic mate- rial. In treating larger bone defects, the material can be combined with autogenous or allogenic bone tissue; the use of membranes depends on the type of operation carried out, the indications, the examining doctor and the procedure recommended in the package leaflet of the given membrane. Consistent oral hygiene both prior to and after dental surgery is necessary for achieving good results. The type of local anaesthetic should be decided by the treating doctor and depends on the package leaflet of the manufacturer.

Antibiotic shield: An antibiotic shield should be used for a period of 7 days. Healing without an antibiotic shield can be tolerated when the correct operation technique is employed and small defects of up to 0.5–1.0 cm are filled in. Dehiscence of the operation wound is an indication for prolonging the antibiotic shield of the operation up to the re-epitheliazation of the defect.

Size of the bioceramic particles: The optimum size of the

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particles is 1/10 of the diameter of the defect. A large bone defect filled with fine ceramic granules can lead to a nonvascularized X-ray contrast block with very low oxygen perfusion that is metabolically completely inactive.

Packaging: The granules are supplied sterile in transparent, glass bottles, closed by a silicon stopper and flip-off seal and put into the blister with the paper cover. The product is supplied with the instruc- tions for use. The product PORESORB-TCP is supplied sterile and is intended for single use only. Reuse of single-use devices creates a risk of patient or user infection. LASAK does not accept any responsi- bility for the re-sterilization of the medical device, regardless of who has carried out re-sterilization or by which method. Do not use if the package is damaged. Non-sterile product may, under no circum- stances, be used. Product must never be cleaned or re-sterilized and must be disposed of. Disposal is not subject to special requirements. The product may not be used after the expiry date.

Storage: The product should be stored in the original packaging, at room temperature and in a dry place.

Copyright: LASAK documentation may not be reproduced by any means without the written permission of LASAK. PORESORB-TCP is a registered trademark of LASAK.

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