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Characteristics: An inorganic nano and macro porous bone regeneration material with osteoconductive properties intended for replacement of bone tissue. The material is white in color with a density of 2,750–3,140 kg/m³. The material bonds directly to bone without any intermediate fibrous tissue layer.

Indications: Replacement of missing or lost bone tissue after insertion into the place of the defect independently or in combination with live bone tissue. The material can be used to shape the bone profile in membrane-covered defects. In dentoalveolar surgery for “sinus lift” operations, augmentation of the atrophic alveolar ridge, filling of bone defects following extirpation of cysts or surgical extractions and after extraction for the prevention of alveolar ridge atrophy. In periodontology, for treatment of one- or multi-walled bone periodontal defects. In implantology, for substitution of the bone defect perioperatively or in peri-implantitis treatment. In orthopedics and traumatology, for treatment of tumor-like affections (bone cysts), bone gangliomas, fibrous dysplasia, pathological fractures, post-traumatic bone defects, benign bone tumors, and arthrodesis. In spinal surgery, for filling bone defects and intervertebral cages.

Contraindications:

a) **Local:** Acute or chronic inflammatory process in the operated area, radiotherapy in the operated area, pathological changes in the jaws or buccal mucosa in dental indications.

b) **Temporary:** Acute febrile disease, pregnancy, treatment with corticosteroids or other drugs affecting calcium metabolism, cytostatic or anticoagulant treatment, extreme physical or mental stress.

c) **General medicinal:** Endocrinopathy and other diseases affecting calcium metabolism, uncontrolled diabetes mellitus, hematological diseases and bleeding disorders, systemic bone diseases, systemic rheumatic diseases, congenital and acquired immune disorders, cytostatic, immunosuppressive or anticoagulant treatment.

Undesirable effects: Not observed to date.

Interactions: None observed to date.

Means of application: Prior to the application of OssaBase-HA, the granules may be mixed with the patient's own blood from the operation wound or withdrawn venous blood, with the patient's own platelet rich plasma (PRP). Application and final shaping of the mixture at the site of the defect is possible using a suitable surgical instrument. The adapted and condensed ceramic is enclosed by the soft tissues using a hermetic suture. It is not advisable to saturate the granules with aqueous solutions (such as normal saline).

Operation technique: OssaBase-HA should be applied by, or under the supervision of, a physician with surgical and/or periorodontal training and experience in the use of biomaterials. The use of the OssaBase-HA material requires gentle operation techniques and prevention of contamination of the material and the operation wound (for example, by saliva in dentoalveolar surgery). Undis-

turbed healing of the OssaBase-HA material requires close contact with the vital bone tissue and assurance of stability in the bone socket. Micro movements of the solid particles cause proliferation of connective tissue between the solid particles and their fibrous encapsulation. The opening of the soft tissues should be extensive, enabling good access to the operation field. It is necessary to uncover the vital bleeding bone tissue for subsequent healing with the ceramic material. In treating larger bone defects, the material can be combined with autogenous or allogenic bone tissue; the use of membranes depends on the type of operation carried out, the indications, the examining doctor, and the procedure recommended in the package leaflet of the given membrane. Consistent oral hygiene, both prior to and after dental surgery, is necessary for achieving good results. The type of local anesthetic should be decided by the clinical doctor and depends on the package leaflet of the manufacturer.

Antibiotic shield: An antibiotic shield should be used for a period of 7 days. Healing without an antibiotic shield can be tolerated when the correct operation technique is employed and small defects of up to 0.5–1 cm are filled in. Dehiscence of the operation wound is an indication for prolonging the antibiotic shield of the operation up to the reepithelialization of the defect.

Size of the bioceramic particles: The optimum size of the particles is 1/10 of the diameter of the defect. A large bone defect

filled with fine ceramic granules can lead to a nonvascularized X-ray contrast block with very low oxygen perfusion that is metabolically completely inactive.

Packaging: The granules are supplied sterile in transparent, sterile bottles, closed by a silicon stopper and flip-off seal. The product is placed into a plastic blister with a paper cover or into the paper box with the instructions for use. The product, OssaBase-HA, is supplied sterile and is intended for single use only. Reuse of single-use devices creates a risk of patient and/or user infection. LASAK does not accept any responsibility for the re-sterilization of the medical device, regardless of who has carried out re-sterilization, or by which method. Do not use if the package is damaged. Non-sterile product may, under no circumstances, be used. Product must never be cleaned or re-sterilized and must be disposed of. Disposal is not subject to special requirements. The product may not be used after the expiry date.

Storage: The product should be stored in the original packaging, at room temperature and in a dry place.

Copyright: LASAK documentation may not be reproduced by any means without the written permission of LASAK. OssaBase-HA is a registered trademark of LASAK.

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